



Product Catalogue

2026/27

Dedicated Team - Free Demos - Training Support - Expertise on Hand - Nationwide Support

Behind the brand

Here at FormAlign our team has its roots right at the heart of the healthcare industry. For over 30 years we have worked closely with healthcare professionals, providing care equipment and have developed a strong reputation to provide the right advice and the right product solution for the service user, based on a depth of understanding of the needs that exist.

Today we focus on postural seating products and component solutions, supporting our network of distributors, wheelchair services and community loan stores throughout the UK. Through close experiences and extensive networking, we have developed products focused around the needs of the industry. Our range of postural components and solutions for wheelchair seating, and cushions form a collection which delivers genuine value.

We know how important it is to deliver the right solutions on time, every time. Our distribution and warehousing operation is set up to offer competitive delivery times to support our trade partners – ensuring the right solution is prescribed on time for the end user who is at the centre of all that we do.

Read our customer Testimonials



My client has said it is the most comfortable and supported she has ever felt sitting in her wheelchair.

Wheelchair Prescriber



WOW its like sitting on a cloud!

End user



Your new gear is impressive, love it!

Wheelchair Prescriber



It's an absolute dream - it is so comfortable - I was able to go for an 8 mile trip recently in my wheelchair and the pain I had in my previous seating, has gone! I am also sleeping better at night - thank you

End user



The combination of the cushion and lateral support has worked exceptionally well to correct the right hip abduction. We have observed significant improvements in both, alignment and comfort, which have greatly benefited our service user. The effectiveness of your products in achieving these outcomes has exceeded our expectations.

Wheelchair Prescriber



Having the support from the FormAlign guys present at the clinic, ensured the best outcome for our service user. And it also gave me a chance to learn another trick or two for their seating.

Wheelchair Prescriber

Contents

Home & School 04

→ B50 Hi-Lo Activity Chair

Cushions 06

→ Pressure Map	07
→ Assist Cushion	08
→ Contour Cushion	08
→ F1 Cushions	09
→ F2 Cushions	09
→ F2X Cushions	09
→ F3 Cushions	09
→ F4 Cushions	10
→ F5 Cushions	10
→ F5+ with Gel Cushions	10
→ STX Cushion	11
→ Woven Air	12
→ Woven Air Base	13
→ Sprint Hybrid Cushion	13
→ Air Insert	14
→ STX Air Cushion	14
→ Air Tubes	14
→ S5 Cushion	15
→ G7 Modular Cushion	15
→ Cushions Rigidiser	16
→ Cushions Lat. Support	16
→ Seat Pan	16

Backrests 17

→ Ignite Backrest	18
→ The Dynamic Bracket Kit	20
→ Contour Backrest	21
→ Support Backrest	21
→ FormAlign Growth Kit	22
→ U2 Second Base Clamps	22
→ 4 Point Upgrade Kit	22
→ U2 Replacement & Mounting Kit	22

Postural Support 23

→ Teardrop Lat Pad	24
→ Straight Lat Pad	24
→ Curved Lat Pad	24
→ Clamp-On Lat Bracket	25
→ Interface Kits	25
→ Height Adjustable Clamp-On Lat	25
→ Fixed Lat Bracket	26
→ Standard Swing-Away Bracket	26
→ Lo-Rise Swing- Away Bracket	26
→ Lo-Rise Swing-Away Bracket	26
→ Axis Headrest Stem Hardware	27
→ Axis Headrest w/ Contour Pad	27
→ Axis Headrest w/Pivot-Wing Pad	27
→ H1 Headrest	28
→ Trip Saver Kit	28

Contents

Positioning and Accessories 29

→ P8 Self-Adhesive Gel Sheet	29
→ P7 Gel Sheet	30
→ FormaAlign Gel Wrap	30
→ Positioning Cluster Pad	30
→ BAC System	31
→ Thigh/Hip Guides	31
→ Swing Down Pommel	31
→ P12 Positioning Kit	32
→ Obliquity Pads	32
→ Self Adhesive Hook & Loop	32
→ Crumb Foam Footplate Pillow	33
→ T-roll	33
→ FormAlign P6	33
→ FormAlign P1	34
→ FormAlign P2	34
→ FormAlign P3	34
→ FormAlign P5	34

Belts & Harnesses 35

→ 2 Point Pelvic Belt	36
→ 4 Point Pelvic Belt	36
→ Chest Harness with Zip	36
→ Chest Harness without Zip	36
→ Shoulder Harness	36
→ Groin Harness	36
→ Calf Strap	36
→ Ankle Huggers	36
→ D-rings	36
→ Universal Strap	36
→ Hardware Tube Clamp	37
→ Cam Buckles	37
→ Tri Glides	37
→ Harness Retainer Kit	37

Home & School Seating

ONE BASE
MULTIPLE CONFIGURATIONS
LONG TERM VALUE

NEW

B50 Hi-Lo Activity Chair

Many children achieve their best postural control in their wheelchair seating — yet once they move into home or school environments, that support often disappears.

The result?

Compromised posture, reduced engagement, and more time spent readjusting than participating.

The HiLo Activity Chair Base was created to extend consistent, clinically effective seating beyond the wheelchair, into the places where learning, play and daily activities actually happen.



Same seating - Same postural outcomes - More environments

COMPONENT BASED



HiLo Base
FGA3001



**Split footplates with
mounting hardware**
FGA3009



Activity Tray
FGA3003



**Armrest LHS Complete
with hardware**
FGA3005



**Backrest and push
handle**
FGA3007



**Central Footplate with
mounting hardware**
FGA3002



Device tray
FGA3008



**Armrest LHS Complete
with hardware**
FGA3006

- Specifically engineered for paediatric needs to provide comfort and stability.
- A design that enables effortless fitting and the ability to re-issue components for long-term sustainability.
- Every feature is designed to be easily adjusted, making the seating experience seamless for the carer.



Cushions

Gel, air, or foam—which one is actually right for your client?

Choosing a cushion doesn't have to be complicated. FormAlign® has removed the guesswork by carefully selecting high-quality materials and combining them into a streamlined range of modular cushions. The result is a more consistent prescribing experience, greater confidence for clinicians, and improved comfort and dignity for users.


Our cushions incorporate the proven benefits of three core technologies—Gel, Air, and Foam—each offering unique clinical advantages. Rather than asking clinicians to choose between them, we've engineered modular solutions that combine the best characteristics of each!

Designed to adapt to a wide range of users and changing clinical needs, our cushions are maintenance-free and require minimal adjustment, saving valuable clinical time while delivering reliable performance


Bubble Foam Technology

Pressure Map Study

When independently tested at a rehabilitation institution, the FormAlign® cushion displayed far better pressure distribution than traditional contoured pressure relieving foam cushions. It exhibited minimal variance when alternating the pelvic tilt angle between posterior and anterior. When mapped the cushion was more comparable to an air cushion than its foam counterparts.



FormAlign® Cushion



"Other" Pressure Relieving Foam Contoured Cushion


Shear Protection

The FormAlign® cushion's "bubble contouring" curved top surface gently restricts the forward tilt of the pelvis, which naturally reduces the risk of shear pressure regardless of the sitting position.



Improved Circulation

The mild forward curvature of the cushion relieves constriction to blood vessels behind the knee helping to improve circulation to the lower limbs.



Positioning For Asymmetrical Postures

The detachable pads shown below are used for accommodating or correcting asymmetrical postures. The pads may be layered according to the degree of the asymmetry and can be accurately positioned anywhere on the base of the cushion. You can cut these to shape with a sharp knife or scissors.



Assist Cushion

The Assist cushion is a comfortable seating solution for long term seated patients that have good core stability. This cushion contours around your patient for long lasting comfort whilst maintaining enhanced weight distribution and improves circulation in their lower limbs.

- ✓ 4 Way Stretch water-resistant cover
- ✓ Includes a set of 4 positioning pads for pelvic obliquity
- ✓ 15mm high contouring
- ✓ Pressure relief
- ✓ 136kgs SWL

Seat Width (")	Seat Depth (")	Code
8	8	FGA1306
	10	FGA1307
10	10	FGA1308
	12	FGA1309
12	12	FGA1310
	14	FGA1311
14	14	FGA1312
	16	FGA1313
	18	FGA1314
16	16	FGA1315
	18	FGA1316
	20	FGA1317
18	16	FGA1318
	18	FGA1301
	20	FGA1319
20	18	FGA1320
	20	FGA1302



Contour Cushion

The Contour cushion offers excellent weight distribution and stability. It is comfortable with its deep contour which assists with hip and thigh alignment. This is the perfect solution for long term seated patients who have reduced core stability and skin integrity.

- ✓ 4 Way Stretch water-resistant cover
- ✓ Includes a set of 4 positioning pads for pelvic obliquity
- ✓ 60mm high contouring*
- ✓ Pressure relief
- ✓ 136kgs SWL

Hip Width (")	Seat Width (") *external cushion width measurement	Seat Depth (")	Code
7	8	8	FGA1112
		10	FGA1077
9	10	10	FGA1113
		12	FGA1078
11	12	12	FGA1062
		14	FGA1066
13	14	14	FGA1067
		16	FGA1068
		18	FGA1114
15	16	16	FGA1069
		18	FGA1070
		20	FGA1071
17	18	16	FGA1072
		18	FGA1073
		20	FGA1074
19	20	18	FGA1075
		20	FGA1076



**Talk to us about the 100mm Contour*



F1 Cushions

Low profile, the F1 pressure care cushion is ideal for community applications. The gel component allows for pressure care to be added to standard riser recliner chairs and domestic seating in the community without compromising the seat size and dynamics.

- ✓ Size: 17 x 17 x 1"
- Code: FG0001

This cushion comes with a 4 way stretch wipeable cover

F2 Cushions

The F2 static cushion is created from a single piece foam and is the entry level foam cushion on the series. Complete with a PU 4 way stretch cover, 2 sided zip and large zip cover to assist with infection control, this simple cushion provides a stable seat for those with limited mobility.

- ✓ Paediatric and adult sizes available
- ✓ Size: 17 x 17 x 4"
- ✓ Plus size options available
- Code: FG0010

This cushion comes with a 4 way stretch wipeable cover




F2X Cushions

Designed for superior comfort through greater immersion, the air flow design is a first within a castellated cushion family, which also assists with temperature control. The 4 way air flow allows effective control of body temperature, combined with the breathable cover. The 4 way stretch cover also means pressure care is high on the agenda.

- ✓ Paediatric and adult sizes available
- ✓ Plus size options available

Size	Code
15x15"	FG0020
16x16"	FG0021
16x18"	FG0022
17x17"	FG0023
18x18"	FG0024
18x20"	FG0025
20x18"	FG0026

This cushion comes with a 4 way stretch wipeable cover



F3 Cushions

The added layer of column gel provides superior seating for stability and pressure care. The unique layered design provides an effective solution for high risk pressure needs. The PU 4 way stretch cover and zip cover also assists with infection control. Particularly effective when temperature control of the user is a challenge and the gel provides shear reduction to maintain increased levels of posture.

- ✓ Paediatric and adult sizes available
- ✓ Plus size options available

This cushion comes with a 4 way stretch wipeable cover



Size	Code
15x15"	FG0030
16x16"	FG0031
16x18"	FG0032
17x17"	FG0033
18x18"	FG0034
18x20"	FG0035
20x18"	FG0036





F4 Cushions

The visco topped cushion in the F4 moulds to the shape of the user as the body warms the visco to allow for an increased positive immersion. The increased contact surface in turn reduces the risk to any specific area of skin by spreading the body loading. This medium level pressure cushion is an effective solution for those at risk of pressure sores.

- ✓ Paediatric and adult sizes available
- ✓ Plus size options available

This cushion comes with a 4 way stretch wipeable cover

Size	Code
15x15"	FG0040
16x16"	FG0041
16x18"	FG0042
17x17"	FG0043
18x18"	FG0044
18x20"	FG0045
20x18"	FG0046

F5 Cushions

The F5 offers low to medium level pressure care whilst allowing effective guidance for thigh alignment. The gentle postural control in the F5 cushion provides support without the restriction access, function or comfort for the semi- active user.

- ✓ Paediatric and adult sizes available
- ✓ Plus size options available

This cushion comes with a 4 way stretch wipeable cover

Size	Code
15x15"	FG0050
16x16"	FG0051
16x18"	FG0052
17x17"	FG0053
18x18"	FG0054
18x20"	FG0055
20x18"	FG0056



F5+ with Gel Cushions

The combined solution of postural control and gel layer works to encourage an effective posture along with a high pressure care solution. The deep thigh and hip guides being integral allow for increased postural alignment. The wipeable cover and 4 way stretch means it is practical for use in many applications.

- ✓ Paediatric and adult sizes available
- ✓ Plus size options available

This cushion comes with a 4 way stretch wipeable cover

Size	Code
15x15"	FG0060
16x16"	FG0061
16x18"	FG0062
17x17"	FG0063
18x18"	FG0064
18x20"	FG0065
20x18"	FG0066



STX Cushion



Looking for structure and adjustability within a cushion?

The STX cushion offers a great deal of structure which in turn makes it one of the most supportive cushions on the market for both symmetrical and severely asymmetrical postures. The modular design allows many asymmetric postural challenges to be accommodated or corrected such as wind-sweeping, pelvic rotation, and range of movement limitations within the hip - all within the one system!

The STX has optional Cluster Gel" packs which can be used to further distribute pressure accurately for your more complex postures and skin integrity challenges. From lowering of femurs to offsetting lateral thigh guides, the options are endless!

- ✓ Core structural design withstands constant forces
- ✓ Modular design – customisable off the shelf solution
- ✓ Integral Ischial well adjustment
- ✓ 4 Way Stretch water-resistant cover
- ✓ 80mm high contouring
- ✓ Able to be reset and re-issued
- ✓ Pressure relief options
- ✓ 136kgs SWL



Hip Width (")	Seat Width (") <i>*external measurement of cushions</i>	Seat Depth (")	Clusters Required	Code
10.5/11	12	12	3	FGA1513
		14	3	FGA1514
12.5/13	14	14	4	FGA1515
		16	4	FGA1516
		18	6	FGA1517
14.5/15	16	16	4	FGA1518
		18	6	FGA1519
		20	6	FGA1520
16.5/17	18	18	6	FGA1521
		20	6	FGA1522
		22	6	FGA1523
18.5/19	20	20	6	FGA1524

17x17" (clusters required 4) FGA1728

Clusters	Code
Foam Cluster	FGA1337
Visco Cluster	FGA1412
Gel Cluster	FGA1238

Cushion sections are available to purchase separately.

Air Cushions

Engineered to reduce pressure and improve ventilation, our air cushions provide lightweight, responsive support for enhanced comfort and long-term skin protection.

Woven Air


The well proven Woven Air. Providing superior posture control that you've come to know, combined with the superior levels of pressure care

- ✓ Maintenance free air system
- ✓ Provides very high levels of pressure care
- ✓ Lateral stability and high-level posture care
- ✓ Available as a component-based unit for third party applications
- ✓ SWL: 180kg



Initial

Across 30 mins



Seat Width (")	Seat Depth (")	Code
14	14	FGA1739
14	16	FGA1740
14	18	FGA1741
15	15	FGA1742
15	17	FGA1743
15	19	FGA1744
16	16	FGA1745
16	18	FGA1746
17	17	FGA1747
17	19	FGA1748
18	18	FGA1749
18	20	FGA1750
18	22	FGA1751
19	19	FGA1752
20	20	FGA1753
20	18	FGA1754

Patented Woven Air™ Technology

How does it work?

Air travels through multiple air-tight tubes which are woven through a tube mat. The air tubes compress under the weight of the user redistributing the pressure.



Woven Air Tube

The result is an adjustable air cushion which has fantastic stability, great comfort, immersion & pressure redistribution without any ongoing maintenance.

The lightweight design offers a great choice for active users. Typically a 18x18 Woven Air cushions weighs in at approx. 1.35kg

Woven Air Inserts are now available.

**Woven Air Tubes are available to purchase separately*



Making Your Day!

WA Positioning Base

Through decades of experience we know that the principles of seating don't change. Compromising on anything is not part of the equation when it comes to FormAlign.

The positioning base for the Woven Air is a useful tool when looking at specific positioning which can be added to the existing Woven Air cushion without the need for a specific cover to suit the additional part.

Whether exploring solutions for a pelvic obliquity, providing additional height in the cushion or accommodating limited hip range of movement, this product delivers options for the user.



Seat Width (")	Seat Depth (")	Code
14	14	FGA1780
14	16	FGA1781
14	18	FGA1782
15	15	FGA1783
15	17	FGA1784
15	19	FGA1785
16	16	FGA1786
16	18	FGA1787
17	17	FGA1788
17	19	FGA1789
18	18	FGA1790
18	20	FGA1791
18	22	FGA1792
19	19	FGA1793
20	20	FGA1794
20	18	FGA1795

- ✔ **Pressure Relief** – Doesn't compromise immersion into air tubes.
- ✔ **Stability** - Firm base provides secure seating.
- ✔ **Custom options** – Can be customised in house or in clinic with the patient. **Easy Transfers** - Mild lateral contours for smooth movements.
- ✔ **Durable** - Maintains support with regular use.

Sprint Hybrid Cushion

Hybrid foam/air

Made for active users!

This ultra-lightweight hybrid cushion has an innovative air/foam hybrid design for maximum stability and energy transfer. Air tubes in the ischial region prevent pressure injuries, and its slim profile lowers the centre of gravity. Plus, it offers top- notch performance with zero maintenance.



Seat Width (")	Seat Depth (")	Code
15	15	FGA2085
	17	FGA2086
	19	FGA2087
16	16	FGA2088
	18	FGA2089
	20	FGA2090
17	17	FGA2091
18	18	FGA2092
	20	FGA2093
	22	FGA2094
20	20	FGA2095

- ✔ **Ultra-Lightweight**: Ideal for active users.
- ✔ **Air/Foam Hybrid**: Combines air and foam for stability and comfort.
- ✔ **Pressure Protection**: Air tubes help reduce pressure injury risk.
- ✔ **Maintenance-Free**: Sealed air tubes require no inflation.
- ✔ **Low Profile**: 55 mm design promotes a lower centre of gravity.
- ✔ **Enhanced Stability**: Foam front helps reduce pelvic forward migration.
- ✔ **Easy Transfers**: Low side contours improve access.
- ✔ **Customisable**: Compatible with positioning pads for additional support.

Air Insert

Our all new Air Insert product is based on our very-high-risk Woven Air cushion which is well proven and providing the market with unparalleled results for prescribers and users alike.

Our insert is now on offer to the market for third party applications too such as:

- ✓ Custom moulded seating systems
- ✓ Foot protection via a wheelchair footplate
- ✓ Backrest options
- ✓ Armrest options

True to our portfolio, FormAlign's truly component based range is designed for the prescriber to create and deliver the outcome required.

Code: FGA1732

Dimensions		
Width	Dept	Height
12"	8.5"	2"



STX 2 Air Cushion

STX Air builds on the proven structure and stability of the original STX cushion, now enhanced with an optional air insert for superior pressure care. Designed for users at very high risk of pressure injury, its advanced air system provides exceptional pressure redistribution while maintaining excellent posture control and stability.

This modular cushion can be finely adjusted over time to meet changing clinical needs, helping users who may have previously been unable to sit out for extended periods improve their comfort, independence, and quality of life.

- ✓ Integrated air insert for high-level pressure care.
- ✓ High lateral and adductor pommel support for pelvis, hips, and thighs.
- ✓ Fully adjustable contours to support symmetrical and asymmetrical sitting positions.
- ✓ Fluid gel technology to redistribute weight and relieve pressure on at-risk areas.
- ✓ Comfort wrap for resilient comfort and cluster gel fluid packs for pressure equalisation.
- ✓ Adjustable positioning pads to address pelvic obliquity, tilt, rotation, femoral rotation, and wind-sweeping.
- ✓ Removable Thigh Wedges: Included to address hamstring shortening and other specific needs.
- ✓ Firm unique one-piece-base base for stability for maintenance of position
- ✓ Built for the prescriber, allowing for component-based adjustments and long-term solutions.

Seat Width (")	Seat Depth (")	Code
15	15	FGA2074
	17	FGA2075
	19	FGA2076
16	16	FGA2077
	18	FGA2078
	20	FGA2079
17	17	FGA2080
18	18	FGA2081
	20	FGA2082
17	22	FGA2083
20	20	FGA2084



Air Tubes

Ideal for those requiring advanced pressure and posture care, the STX Air Cushion is the ultimate solution for enhanced comfort and support.



Size	Code
14"	FGA1755
15"	FGA1756
16"	FGA1757
17"	FGA1758
18"	FGA1759
20"	FGA1760
22"	FGA1761

S5 Cushion

For our users with high pressure care needs

Superior in pressure care, the pre-moulded SX5 cushion offers great comfort with its soft, hybrid gel system, making it a great choice for wheelchair users who do not have complex postural needs but providing stability, high level pressure care and ultimate comfort. It's gentle, supportive feel is also perfect for aged care, providing lasting comfort for long-term sitting.

- ✓ **Optimal for Aged Care:** Soft, comfortable design perfect for long-term sitting needs.

This cushion combines innovative design with comfort, making it the perfect choice for users seeking soft, supportive seating without postural asymmetry concerns.

Seat Width (")	Seat Depth (")	Code
15	15	FGA2108
	17	FGA2109
	19	FGA2110
16	16	FGA2111
	18	FGA2112
	20	FGA2113
17	17	FGA2114
	18	FGA2115
	20	FGA2116
18	22	FGA2117
	20	FGA2118
20	20	FGA2118



Soft Foam Base - Adds extra comfort for extended sitting periods.



Fluid "Gel Pax" Technology - Ensures even pressure dispersion and enhanced comfort.



(Gel pads are inserted into the integrated pockets)



Cover Included

G7 Modular Cushion

Created for outstanding stability, comfort, and pressure redistribution, the G7 cushion is ideal for users with asymmetrical postures and at very high risk of pressure sores. Its mild lateral contouring facilitates easy side and frontal transfers, enhancing overall ease of use.

- ✓ **Pressure Relief** - Redistributes pressure away from high-risk areas via gel pax system.
- ✓ **Stability** - Firm base provides foundation for seating.
- ✓ **Customisation** - Adjustable pads included for pelvic alignment.
- ✓ **Easy Transfers** - Mild lateral contours for smooth movements.
- ✓ **Durable** - Maintains support with regular use.
- ✓ **Washing** - Removable, machine-washable cover.


Seat Width (")	Seat Depth (")	Code
15	15	FGA2060
	17	FGA2061
	19	FGA2062
16	16	FGA2063
	18	FGA2064
	20	FGA2065
17	17	FGA2066
	18	FGA2067
	20	FGA2068
18	22	FGA2069
	20	FGA2070
20	20	FGA2070



Cushions Rigidiser

The cushion rigidiser is designed to work in combination with a cushion on a wheelchair seat canvas to prevent "sagging" and loss of structure to ensure best postural control. It is simply placed inside the cushion. The composite material which this product is manufactured from maintains structure whilst keeping the weight to a minimum.

- ✓ 4mm thick
- ✓ Velcro strips for secure adhesion to base



Seat Width (")	Seat Depth (")	Code
8	8	FGA1157
	10	FGA1158
10	10	FGA1159
	12	FGA1160
12	12	FGA1129
	14	FGA1130
14	14	FGA1131
	16	FGA1132
	18	FGA1133
16	16	FGA1134
	18	FGA1135
	20	FGA1136
18	16	FGA1137
	18	FGA1138
	20	FGA1139
20	18	FGA1161
	20	FGA1162

Cushions Lateral Support

This innovative yet simple product has been designed to provide extra support when required against strong abduction.

- ✓ Fits inside cushion cover
- ✓ Velcro strips for secure fit
- ✓ Lightweight



Lateral Support Height (")	Seat Width (")	Code
3	8	FGA1115
3	10	FGA1116
3	12	FGA1117
3	14	FGA1118
4	16	FGA1080
4	18	FGA1119
4	20	FGA1120

Measurements:

Height of lateral: 3"

Seat Pan

KEPT IN STOCK

The quick release seat pan is a great solution for use on folding mobility chassis. Transport safety approved to WC20 and ISO16840-4 standards for forward facing transportation in a motor vehicle.

The seating system can be mounted onto the seat pan and removed easily for storage by parents and carers.

- ✓ Ideal for wheelchair storage restrictions
- ✓ Fully crash tested



To fit to rails 17-25mm diameter

Seat Width (")	Seat Depth (")	Code
12	12	FGA1104
	14	FGA1105
14	14	FGA1106
	16	FGA1081
	18	FGA1107
16	16	FGA1108
	18	FGA1055
	20	FGA1109
18	18	FGA1110
	20	FGA1054
20	20	FGA1111

Backrests

A young boy with short brown hair is seen from the side, sitting in a wheelchair. He is wearing a light blue t-shirt and dark shorts. He is looking out over a grassy field towards some trees in the distance. The background is slightly blurred. A large orange circle is overlaid on the left side of the image, containing the word 'Backrests' in white text.

Ignite



Can be used anyway up to achieve multiple thoracic lateral height position

Ignite Backrest

The Ignite Backrest is designed for the active user and single shell backrest market. With a aluminium shell as well as the proven "bubble foam" for breathability & comfort as well as the U2 hardware, we have you backed all the way!

- ✓ Available in sizes from Paediatric to Adult
- ✓ Easy to prescribe, easy to adjust and easy to order!
- ✓ Can be recycled and reissued

Comes complete with mounting hardware. SWL: 136kg

DEEP Contour Shoulder Height Backrest complete with quick release hardware

User Chest Width (")	User Back Height (")	Backrest Size (")	Will fit Wheelchair Width (")	Lateral Support Depth (")	Will fit Backrest Tube Dia. (")	Code
11	16	14 x 16	12 up to 16	5	Up to 1	FGA1642
11	18	14 x 18	12 up to 16	5	Up to 1	FGA1643
12	17	15 x 16	13 up to 17	5	Up to 1	FGA1644
12	17	15 x 18	13 up to 17	5	Up to 1	FGA1645
13	18	16 x 18	14 up to 18	6	Up to 1	FGA1646
13	20	16 x 20	14 up to 18	6	Up to 1	FGA1647
14	18	17 x 18	15 up to 19	6	Up to 1	FGA1648
14	20	17 x 20	15 up to 19	6	Up to 1	FGA1649
15	20	18 x 18	16 up to 20	6	Up to 1	FGA1650
15	20	18 x 20	16 up to 20	6	Up to 1	FGA1651
17	18	20 x 18	18 up to 22	6	Up to 1	FGA1652
17	20	20 x 20	18 up to 22	6	Up to 1	FGA1661



From the comfortable pressure re-distributing foam inner to the easy to use U2 compact quick-release hardware, the light weight Dreamline Ignite Backrest System is also available in three contour depths and three different heights. FormAlign stocks the **Deep Contour**, **Shoulder Height** backrest but other depths and heights are available to order on request.

CONTOUR DEPTHS



DEEP CONTOUR



MID CONTOUR



ACTIVE CONTOUR

SHELL HEIGHT OPTIONS



SHOULDER HEIGHT



THORACIC HEIGHT



LUMBAR HEIGHT

CONTOUR TYPE	HEIGHT	USER CHEST WIDTH (")	LATERAL DEPTH	BACKREST SHELL HEIGHT (")	BACKREST SIZE WIDTH x HEIGHT (")	Wheelchair Frame Width (Outside to Outside)	Physical Backrest Shell Width	Hardware Clamps (Tube Diameter)	ORDER CODE
ACTIVE	Lumbar	9-10"	1.5"	7"	14" x 7"	12" to 16"	11"	Up to 1"	FGA1686
ACTIVE	Thoracic	9-10"	1.5"	12"	14" x 12"	12" to 16"	11"	Up to 1"	FGA1687
ACTIVE	Shoulder	9-10"	1.5"	16"	14" x 16"	12" to 16"	11"	Up to 1"	FGA1688
ACTIVE	Lumbar	10"-11"	1.5"	7"	15" x 7"	13" to 17"	12"	Up to 1"	FGA1689
ACTIVE	Thoracic	10"-11"	1.5"	12"	15" x 12"	13" to 17"	12"	Up to 1"	FGA1690
ACTIVE	Shoulder	10"-11"	1.5"	16"	15" x 16"	13" to 17"	12"	Up to 1"	FGA1691
ACTIVE	Lumbar	11"-12"	2"	8"	16" x 8"	14" to 18"	13"	Up to 1"	FGA1692
ACTIVE	Thoracic	11"-12"	2"	13"	16" x 13"	14" to 18"	13"	Up to 1"	FGA1693
ACTIVE	Shoulder	11"-12"	2"	18"	16" x 18"	14" to 18"	13"	Up to 1"	FGA1694
ACTIVE	Lumbar	12"-13"	2"	8"	17" x 8"	15" to 19"	14"	Up to 1"	FGA1695
ACTIVE	Thoracic	12"-13"	2"	13"	17" x 13"	15" to 19"	14"	Up to 1"	FGA1696
ACTIVE	Shoulder	12"-13"	2"	18"	17" x 18"	15" to 19"	14"	Up to 1"	FGA1697
ACTIVE	Lumbar	13"-14"	2"	8"	18" x 8"	16" to 20"	15"	Up to 1"	FGA1698
ACTIVE	Thoracic	13"-14"	2"	13"	18" x 13"	16" to 20"	15"	Up to 1"	FGA1699
ACTIVE	Shoulder	13"-14"	2"	18"	18" x 18"	16" to 20"	15"	Up to 1"	FGA1700
ACTIVE	Lumbar	14"-15"	2"	8"	20" x 8"	18" to 22"	16"	Up to 1"	FGA1701
ACTIVE	Thoracic	14"-15"	2"	13"	20" x 13"	18" to 22"	16"	Up to 1"	FGA1702
ACTIVE	Shoulder	14"-15"	2"	18"	20" x 18"	18" to 22"	16"	Up to 1"	FGA1703
MID	Lumbar	9-10"	3"	7"	14" x 7"	12" to 16"	11"	Up to 1"	FGA1721
MID	Thoracic	9-10"	3"	12"	14" x 12"	12" to 16"	11"	Up to 1"	FGA1704
MID	Shoulder	9-10"	3"	16"	14" x 16"	12" to 16"	11"	Up to 1"	FGA1705
MID	Lumbar	10"-11"	3"	7"	15" x 7"	13" to 17"	12"	Up to 1"	FGA1706
MID	Thoracic	10"-11"	3"	12"	15" x 12"	13" to 17"	12"	Up to 1"	FGA1707
MID	Shoulder	10"-11"	4"	16"	15" x 16"	13" to 17"	12"	Up to 1"	FGA1708
MID	Lumbar	11"-12"	4"	8"	16" x 8"	14" to 18"	13"	Up to 1"	FGA1709
MID	Thoracic	11"-12"	4"	13"	16" x 13"	14" to 18"	13"	Up to 1"	FGA1710
MID	Shoulder	11"-12"	4"	18"	16" x 18"	14" to 18"	13"	Up to 1"	FGA1711
MID	Lumbar	12"-13"	4"	8"	17" x 8"	15" to 19"	14"	Up to 1"	FGA1712
MID	Thoracic	12"-13"	4"	13"	17" x 13"	15" to 19"	14"	Up to 1"	FGA1713
MID	Shoulder	12"-13"	4"	18"	17" x 18"	15" to 19"	14"	Up to 1"	FGA1714
MID	Lumbar	13"-14"	4"	8"	18" x 8"	16" to 20"	15"	Up to 1"	FGA1715
MID	Thoracic	13"-14"	4"	13"	18" x 13"	16" to 20"	15"	Up to 1"	FGA1716
MID	Shoulder	13"-14"	4"	18"	18" x 18"	16" to 20"	15"	Up to 1"	FGA1717
MID	Lumbar	14"-15"	4"	8"	20" x 8"	18" to 22"	16"	Up to 1"	FGA1718
MID	Thoracic	14"-15"	4"	13"	20" x 13"	18" to 22"	16"	Up to 1"	FGA1719
MID	Shoulder	14"-15"	4"	18"	20" x 18"	18" to 22"	16"	Up to 1"	FGA1720
DEEP		11"	5"	11"	16x13"	12 to 16"	14"	Up to 1"	FGA1655
DEEP		12"	5"	11"	17x13"	13 to 17"	14"	Up to 1"	FGA1656
DEEP		13"	6"	13"	16x13"	14 to 18"	16"	Up to 1"	FGA1657
DEEP		15"	6"	13"	17x13"	15 to 19"	17"	Up to 1"	FGA1658
DEEP		15"	6"	13"	18x13"	16 to 20"	18"	Up to 1"	FGA1659
DEEP		17"	6"	13"	20x13"	18 to 22"	20"	Up to 1"	FGA1660

The Dynamic Bracket Kit

The Dynamic bracket kit is designed to be added to existing backrests as clinical needs present or change or re-issuing is required. This clever, component-based design allows the seating to be built to meet the users specific needs.

- ✓ No modification to the wheelchair frame
- ✓ Protects the wheelchair equipment from rocking force
- ✓ Uses the same, trusted hardware
- ✓ Adjustable resistance pads
- ✓ Retro-fits to existing bracket shells
- ✓ Minimal maintenance

This product comes with soft, medium and firm elastomer pads which can be easily changed (all included with part number FGA1922)

Why the Dynamic Bracket?

- ✓ 2" of suspension travel.
- ✓ Lock for travel in a motor vehicle or tilt-in-space.
- ✓ Easy access points to quickly change or replace elastomer pads.
- ✓ A height adjustable top mount enables further adjustments to movement resistance.
- ✓ Quick-release 4 point attachment

Code: FGA1922



**Backrest not included in the kit*

Contour Backrest

Contour Backrest is a one-piece shell with integrated lateral supports. The Contour's signature "bubble" improves immersion, enabling the backrest to closely match the contour of the user.

This is perfect for the semi active user, providing comfortable torso support and stability which, in turn, encourages a mid-line sitting position.

- ✓ Can be recycled and reissued
- ✓ Fully Crash Tested
- ✓ Pressure relief
- ✓ Compatible with third party wheelchairs
- ✓ 136kg sSWL

Comes complete with mounting hardware.

"These backrests will fit to wheelchairs that are 2" bigger or 2" smaller than the specified backrest size"



Wheelchair Width (")	Backrest Height (")	Trunk Width (")	Code
8	10	6	FGA1084
	12	6	FGA1088
10	12	8	FGA1085
	14	8	FGA1089
12	14	10	FGA1001
	16	10	FGA1002
	18	10	FGA1003
14	16	12	FGA1004
	18	12	FGA1005
	20	12	FGA1006
16	18	14	FGA1007
	20	14	FGA1008
	22	14	FGA1009
18	18	16	FGA1010
	20	16	FGA1011
	22	16	FGA1012
20	18	18	FGA1013
	20	18	FGA1014
	22	18	FGA1015

Support Backrest

Support Backrest is a perfect platform for attaching modular supports such as swing-away or fixed adjustable lateral supports.

This is a great solution for less active users with symmetrical or asymmetrical postures.

NEW Updates to the Support backrest! As standard we now have the Flo2 Cover which is multi-direction stretch which allows market leading immersion properties and thus superior comfort and stability for the user. Its slim foam inner gives you more seating depth making it a perfect option for asymmetric body positions.

The BAC positioning kit inserts even more subtly when needed and can be added and removed as users needs require and is used for asymmetric postures to promote increased contact and feedback.

The Dreamline Support Backrest combines innovative design with advanced materials to provide unbeatable comfort and support for less active users.

- ✓ Can be recycled and reissued
- ✓ Fully Crash Tested
- ✓ Pressure relief
- ✓ Compatible with third party wheelchairs

Comes complete with mounting hardware.

"BAC system compatible for the Support Backrest. See page 20 for details."

"These backrests will fit to wheelchairs that are 2" bigger or 2" smaller than the specified backrest size"

Wheelchair Width (")	Backrest Height (")	Physical Shell Width (")	Code
8	10	5	FGA1082
	12	5	FGA1090
10	12	7	FGA1083
	14	7	FGA1091
12	14	9	FGA1016
	16	9	FGA1017
	18	9	FGA1018
14	14	11	FGA1163
	16	11	FGA1019
	18	11	FGA1020
16	20	11	FGA1021
	16	13	FGA1164
	18	13	FGA1022
	20	13	FGA1023
	22	13	FGA1024
18	18	15	FGA1025
	20	15	FGA1026
	22	15	FGA1027
20	18	17	FGA1028
	20	17	FGA1029
	22	17	FGA1030

For backrests, under 12", please use Kids Lo Rise Lateral Brackets,



Assymetrical postures





FormAlign® Growth Kit

The FormAlign® Growth Kit enables the backrest to grow in height up to an additional 4".

Backrest Width (")	Code
12	FGA1122
14	FGA1123
16	FGA1124
18	FGA1125
20	FGA1126



U2 Second Base Clamps Hardware

The U2 Second Base Clamps Hardware are for attaching the backrest to a 2nd mobility base.

- ✓ Easy to install and use
- ✓ Crash tested
- ✓ Unique triple lock mount

Clamp size	Code
50mm clamps for Contour/Support Backrest	FGA1167
35mm clamps for Ignite Backrest	FGA1919



4 Point Upgrade Kit

An optional kit to mount the backrest at 4 points as opposed to the standard two points. This can be used where a user may have strong extensor patterns or constant involuntary movements and offers peace of mind from a durability viewpoint.

Code: FGA1140



U2 Replacement & Mounting Kit

These clamps are complete with L brackets.

Kit with Clamp Height	Code
50mm Clamp Height	FGA1154
35mm Clamp Height (Usually used for Ignite systems)	FGA1800

Postural Supports

Support right the way through

All our postural products are component based and have been designed to provide you with the tools to create the solution you need. From building a complete seating system to selecting specific postural parts FormAlign® are here to support you right the way through!

Lateral Pads


Our range of lateral pads can be purchased separately which empowers the prescriber to decide which pad works for each specific case.

Teardrop Lateral Pad

The FormAlign® Teardrop Lateral Pad comes in three different sizing options and can be used on any lateral bracket, shown over.



Pad Size	Dimensions (")	Code
Small	6 H x 6 L	FGA1049
Medium	7.5 H x 6.5 L	FGA1050
Large	9 H x 7.5 L	FGA1051



Straight Lateral Pad

The FormAlign® Straight Lateral Pad comes in three different sizing options and can be used on any lateral bracket, shown over.



Pad Size	Dimensions (")	Code
X-Small	3.5 x 5.5 L	FGA1041
Small	4.5 H x 6.5 L	FGA1042
Medium	5 H x 7 L	FGA1043
Large	5 H x 8.5 L	FGA1044



Curved Lateral Pad

The FormAlign® Curved Lateral Pad comes in three different sizing options and can be used on any lateral bracket, shown over.



Pad Size	Dimensions (")	Code
X-Small	3.5 x 5.5 L	FGA1045
Small	4.5 H x 6.5 L	FGA1046
Medium	5 H x 7 L	FGA1047
Large	5 H x 8.5 L	FGA1048



Cane Mounts

Swing- away brackets are available with our range of lateral pad options to offer the best fit, level of support and stability for the end user. The triple lock mounting hardware secures the clamp without risk of slipping.

FormAlign® Clamp-on Lateral Bracket

The clamp on lateral bracket can be clamped to wheelchair backrest canes/back posts and will fit to cane diameters of 17mm – 25mm (3/4 of an inch to 1") without the need of spacers. The dual pivot bracket delivers market leading adjustment and the support for the user just where it's needed.

- ✔ Will fit to wheelchair when tension adjustable backrest is required
- ✔ Will fit to most brands of wheelchair chassis
- ✔ Durable design
- ✔ Can be recycled and re-issued

Bracket Type	Left/Right	Code
Clamp-on Lateral Bracket	Left	FGA1323
Clamp-on Lateral Bracket	Right	FGA1324



Interface Kits

The interface brackets create solutions to achieve multiple positions for the thoracic support pad for maximum effectiveness. Can be used:

- ✔ On Tension adjustable backrest upholstery
- ✔ When trunk is in forward position due to kyphosis
- ✔ Third party seating frames
- ✔ Restricted arm movement
- ✔ When trying to maximise function such as self-propel motion

Ideal for Self Propel Users



Ideal for Powerchairs



Comes as standard with FGA1323 and FGA1324

Measurements are for a guide only

Plate Size	Dimensions(WxL) (")	Code
Long Plate	3 x 5"	FGA1902
Mid Plate	3 x 2.5"	FGA1901
Short Plate	3 x 1.5"	FGA1900

*When ordering FGA1901 and/or FGA1902, FGA1900 must be ordered also for mounting purposes

Height Adjustable Clamp-On Lateral

Designed for mounting direct to the back posts of the wheelchair, this clamp-on lateral has:

- ✔ Integral height adjustment – no more clashing with hardware
- ✔ Is fully multi-adjustable
- ✔ Will fit posts from 17mm-25mm with no spacers
- ✔ Can be fitted over the top of upholstery
- ✔ Has a pressure relieving gel pad as standard
- ✔ Has hygienic wipeable covers to make it practical and sanitary

Bracket Type	Left/Right	Code
HA Clamp-On Lateral	Left	Available on request
HA Clamp-On Lateral	Right	Available on request

*Please note you can use these brackets as either a left or right



Shell Mounts

Adjustable fixed and swing-away brackets are available with our range of lateral pad options to offer the best fit, level of support and stability for the end user.

FormAlign® Fixed Lateral Bracket

The FormAlign® fixed brackets are width, height and depth adjustable with an off-set attachment position to enable enhanced attachment opportunities.

- ✓ Can be recycled and reissued
- ✓ Compatible with Support Backrest
- ✓ Compatible with third party backrests
- ✓ Designed to meet asymmetric postures

Bracket Type	Code
Fixed Bracket	FGA1469



Standard Swing-Away Bracket

The standard Swing-Away hardware offers a durable system which has an inbuilt dual pivot design, essentially offering an inbuilt "off set" bracket, allowing for maximum trunk support for the user and ultimate positioning. In turn this allows for pressure to be spread over the largest possible surface area which increases comfort by reducing pressure.

- ✓ Can be recycled and reissued
- ✓ Compatible with FormAlign® Support Backrest
- ✓ Compatible with third party backrests
- ✓ Designed to meet asymmetric postures

Bracket Type	Left/Right	Code
Standard Bracket	Left	FGA1052
Standard Bracket	Right	FGA1053



Lo-Rise Swing-Away Bracket

Lo-Rise Swing-Away hardware is the next step up from the standard Swing-Away hardware. The main difference being the vertical adjustment slot. This slot offers superior levels of pad adjustment. It allows us to overcome the well-known challenge of hardware clashes leading to trunk positioning compromises. This solution is very useful when the space is limited on the backrest.

- ✓ Unique integral adjustment and durability
- ✓ Unique dual pivot bracket – negates need for "Z" bracket
- ✓ Twin depth adjustment plate for pad position
- ✓ Compatible with third party backrests

Bracket Type	Left/Right	Code
Lo-Rise Bracket	Left	FGA1039
Lo-Rise Bracket	Right	FGA1040
Lo Rise Lateral Bracket Kids	Left	FGA1141
Lo Rise Lateral Bracket Kids	Right	FGA1142



For backrests, under 12", please use Kids Lo Rise Lateral Brackets

Lo-Rise Swing-Away Pivot Bracket

The additional Pivot variant of the Lo-Rise Swing-Away bracket allows the pad to have further angle adjustment up to 45 degrees-useful for asymmetric trunk shapes and ensuring maximum distribution of load/pressure. This hardware offers the most adjustment and is the most adjustable on the market.

- ✓ Most adjustable in the range with additional vertical Lo-Rise slot
- ✓ Angle adjustment for pad to meet severe scoliosis
- ✓ Unique dual pivot bracket – negates need for "Z" bracket
- ✓ Twin depth adjustment plate for pad position
- ✓ Compatible with third party backrests

Bracket Type	Left/Right	Code
Lo-Rise Pivot Bracket	Left	FGA1384
Lo-Rise Pivot Bracket	Right	FGA1385





Axis Headrest Stem Hardware

The Axis hardware is another word for durability! The design is only a true engineer's handywork. Alongside the heavy-duty system, it is also infinitely adjustable and will offer head support for both kyphotic positions and hyper-extended presentations, along with lateral offset positions to accommodate an asymmetric pattern.

- ✓ 13" forward adjustment as standard
- ✓ Extra links can be added for further reach
- ✓ Lateral positioning
- ✓ Durable mini cog design delivers heavy duty solution
- ✓ Swan neck stem allows for wider range of reach
- ✓ Can be used on third party backrests
- ✓ Complete with mount hardware

Code: FGA1436

Axis Headrest with Contour Pad

The Contour Pad is to support with lateral positioning and can be effective for a guide for the user as well as correction for the head position. Mounted on the Axis stem hardware, this kit is highly adjustable but easy to use daily.

- ✓ Pad available in 2 sizes
- ✓ Occipital shaping for extra support and comfort
- ✓ Foam removal for modifications
- ✓ Self-Adhesive gel available for added pressure relief shear reduction
- ✓ Comes complete with mounting hardware



Pad Size	Dimensions (")	Code
Small	8.5 H x 5.5 W	FGA1127
Large	11 H x 6.5 W	FGA1128



Pad Size	Dimensions (")	Code
Small	5.5 H x 7 W	FGA1303
Large	6.5 H x 6.5 W	FGA1204

Axis Headrest with Pivot-Wing Pad

This kit is an unstoppable solution! The perfect answer for the engineer, the therapist and the repair department – so adjustable, easy to fit and durable! With this one kit, infinite adjustment is an understatement. Paralleled to this, the extremely durable design stands it apart from most headrests for complex needs.

- ✓ Pad available in 2 sizes
 - ✓ Wings are width and depth adjustable
 - ✓ Removable foam sections for pressure care upgrades or environmental controls
 - ✓ Lateral pads can be used any way up for ultimate head positioning
- Self-Adhesive gel available for added pressure relief and shear reduction
Comes complete with mounting hardware



*Pad to be purchased separately.

H1 Headrest

H1 is designed to offer a solution as a guide for the user's head. The contour pad is effective and is used for those with reasonable head control but may require some support during times of fatigue or when using the tilt in space feature.

- ✓ Dual ball pivot system for effective positioning
- ✓ Height and depth adjustment
- ✓ Can be used on third party backrests
- ✓ Complete with mount

Size: Width: 9" Height: 5"

Code: FGA1475



BUILD YOUR OWN TRIPSAVER KIT

Configure a kit that suits your service requirements.

Every clinician works differently. We can configure your TripSaver Kit with the products most relevant to your service, helping you stay organised and prepared wherever assessments take place.

WHY CHOOSE A CUSTOM KIT?



SELECT THE PRODUCTS YOU USE MOST



ORGANISED AND EASY TO TRANSPORT



IDEAL FOR DEMONSTRATIONS AND ASSESSMENTS



DESIGNED AROUND YOUR SERVICE REQUIREMENTS



PRODUCT CODE

FGA1600

FULLY CUSTOMISABLE



Positioning & Accessories

Components to deliver postural support and pressure care solutions

Our handpicked positioning products are all designed to provide the prescriber and rehab dealers the tools to complete the seating system to the best of their ability and for the best user experience and comfort. Whether looking to accommodate an asymmetric foot shape position or to provide added postural stability, our solutions come with many years of hands-on experience and awareness of the daily challenges faced in the specialist seating industry.

P8 Self-Adhesive Gel Sheet

Our FormAlign® Gel sheet product is another clever solution to provide additional pressure care and extra comfort for your customers. Created especially for you, the sheet has a self-adhesive backing and can be cut to the exact shape making it perfect for adding to existing areas of the seating such as headrest pads, trunk support pads and pommel blocks. The sheet is very stretchable and will fold round the contours of an existing postural component which may be in use. The column design allows the surface to reduce shearing on the skin and is effective in maintaining posture. Its extremely lightweight structure enables users to enjoy its benefits whilst hardly even knowing it is there!

Please note product FG9037 is not to be used on top of a seating cushion. Refer to P7 Gel Sheet and Gel Wrap options.

- ✓ Shear reducing design
- ✓ Lightweight
- ✓ Easy to cut and shape
- ✓ Self-Adhesive Backing
- ✓ Non-Absorbing
- ✓ Can be used on endless applications
- ✓ One gel sheet measures 42 x 46cm and is 1.3cm in depth

Code: FG9037



With endless applications, talk to us today about your specific project requirements!


P7 Gel Sheet

This is the solution which offers users the ability to sit longer without discomfort. Lightweight and low profile, the P7 Gel product is designed to work with other seating products such as custom wheelchair seating and contract seating where improved skin care is required. The gel component allows for pressure care to be added into a vast range of seating applications and can be applied to foam seating by seating manufacturers and service providers.

- ✓ High level pressure care
- ✓ Proven pressure redistribution
- ✓ Effective temperature control design
- ✓ Shear reducing design
- ✓ Worry-free comfort

Size: 16x18"

Code: FG9067



FormAlign® Gel Wrap

Our gel wrap system is designed to offer the user extra pressure care. The continuous surface area works well for some active users and those who have specific daily needs. The product can be retrofitted to most cushions and is a simple solution to providing the user with an increased level of pressure care.

Designed to work with the STX cushion, but can be used on other cushions – contact us to discuss your requirements.

- ✓ Range of size options
- ✓ Wipeable outer cover
- ✓ 4 way stretch fabric assists with pressure care

Seat Width (")	Seat Depth (")	Code
12	12	FGA1476
	14	FGA1477
	14	FGA1478
14	16	FGA1479
	18	FGA1480
	16	FGA1481
16	18	FGA1482
	20	FGA1483
	18	FGA1484
18	20	FGA1485
	22	FGA1486
	20	FGA1487



Positioning Cluster Pad

The foam cluster pads are used within the cushion range to create a base within the ischial well for superior immersion. The small clusters are designed to work with adjustable well options and a good choice when working with a custom well and fixed pelvic asymmetry. These clusters can be upgraded to the visco and gel version to increase pressure care levels when required.

- ✓ Can be used on custom ischial wells for accommodating pelvic asymmetries
- ✓ Velcro backing increases application choices
- ✓ 4 way stretch cover
- ✓ Allows for immersion and thus spreads the loading of pressure
- ✓ Can be used in shell backrest products



Pad Type	Code
Foam	FGA1337
Visco	FGA1412
Gel	FGA1238

Each pad measures approximately Size 4x6" in

BAC System

Accommodating mild asymmetrical back profiles is often challenging.

Products that do not offer adjustment to meet these shapes often introduce risks, instability and lack of function for the user. The BAC positioning system has been launched to help to solve these challenges simply with its integral adjustment sections which means mild spinal and rib cage asymmetries can be accommodated easily.

The BAC product is available in a range of sizes and is designed to fit inside the Support backrest. It can be added when the user requires it as needs change and develop.

- ✓ Compatible with most shelled backrests in the market
- ✓ Designed to meet asymmetric postures
- ✓ Provides increased levels of stability and user function
- ✓ Can be reset and re-issued

Size: 18x16"

Code: FGA1410



Backrest Width (")	Backrest Height (")	Physical Width (")	Code
8	10	5	FGA1086
	12	5	FGA1092
10	12	7	FGA1087
	14	7	FGA1093
12	14	9	FGA1064
	16	9	FGA1065
	18	9	FGA1094
14	14	11	FGA1165
	16	11	FGA1063
	18	11	FGA1095
16	20	11	FGA1096
	16	13	FGA1056
	18	13	FGA1057
18	20	13	FGA1058
	22	13	FGA1097
20	18	15	FGA1059
	20	15	FGA1060
	22	15	FGA1098
22	18	17	FGA1099
	20	17	FGA1061
	22	17	FGA1100

Thigh/Hip Guides

The hip guides are easy to adjust in height, depth and width. The hip guides may also be rotated to suit internal/external femoral rotation or wind-swept postures. FormAlign® wheelchair Hip-Guides are transport safety approved and designed to attach to the seat pan enabling precise positioning for the user. The unique flip-down function maintains the pads in a ready position for quick and easy transfers with wide open access to the user.



Code	Dimensions (")
FGA1143	10" x 3" x 1.5"
FGA1178	12" x 3" x 1.5"
FGA1179	14" x 4" x 1.5"
FGA1180	16" x 4" x 1.5"
FGA1144	18" x 4" x 1.5"

- ✓ Adjustable for symmetrical and asymmetrical postures
- ✓ Pads overlaid with comfortable high-density pressure relieving foam
- ✓ Can be rotated to suit internal, external femoral rotation or wind-swept postures.
- ✓ Available in range of sizes


Swing Down Pommel

The swing-down pommel offers comfortable, firm support. With its pressure relieving foam overlay inside a water-resistant cover, it makes a perfect solution for that added thigh support.

- ✓ Compact push button release
- ✓ Height, depth and width adjustable



Pad Size	Dimensions (")	Code
Small	4 L x 3 H x 3 W at front	FGA1147
Medium	6 L x 4 H x 3 W at front	FGA1148
Large	7 L x 4 H x 4 W at front	FGA1149




P12 positioning Kit

Lumbar and lateral wedge kit designed to be purchased for any backrest situation to assist with adjusting the backrest profile for each user. Velcro strips are already adhered to the components making the kit simple to use and configure.

- ✔ Closed-cell foam delivers structure for positioning
- ✔ Velcro strip already adhered
- ✔ Easy to cut to shape
- ✔ Can be used on almost every backrest support

Code: FG1554



Pad Types	Dimensions(")
Side Wedges	4.5 H x 3.5 W
Lumber Pad	6.5 H x 3.5 W

Obliquity Pads

The obliquity detachable pads are designed for accommodating and correcting asymmetrical postures. The pads can be layered according to the degree of the asymmetry, can be accurately positioned anywhere on the base of the cushion and can be cut to shape with a sharp knife or scissors.

Size (") <i>*refer to width of cushion</i>	Dimensions (")	Code
8	4 x 4	FGA1184
10	5 x 5	FGA1185
12	5.5 x 5.5	FGA1186
14	6.5 x 6.5	FGA1187
16	7.5 x 7.5	FGA1188
18	8.5 x 8.5	FGA1189



Self Adhesive Hook & Loop

Our self-adhesive hook and loop solution is again designed with you and your user in mind. The self-adhesive component has an acrylic backing and has an ultra-strong adhesion to many surfaces and is especially effective in the special seating industry.

- ✔ Ultra-strong self-adhesive backing
- ✔ 25mm or 50mm width available
- ✔ Available on 25m roll



Description	Code
50mm wide x 25m roll - LOOP	FGA1491
50mm wide x 25m roll - HOOK	FGA1490

Crumb Foam Footplate Pillow

The Crumb Foam Footplate Pillow is designed to strap onto existing footplates. The Crumb foam filling can be moved to accommodate asymmetric foot shapes. It's the perfect solution for vulnerable skin and high-risk pressure points.

- ✓ Designed to strap onto existing footplates
- ✓ Crumb foam filling can be moved to accommodate asymmetric foot shapes
- ✓ Perfect solution for vulnerable skin and high-risk pressure points

Dimensions(")	Code
15 L x 11 W x 3.5 H	FG9003



T-Roll

Designed to maintain midline leg position when lying supine, the T Roll can be used to prevent lower limb contractures.

- ✓ Designed to maintain midline leg position when lying supine
- ✓ Can be used to prevent lower limb contractures
- ✓ Effective to prevent wind-sweeping
- ✓ Filling can be adjusted to allow for increased immersion
- ✓ Crumb foam filling provides effective pressure care
- ✓ Wipeable covers
- ✓ Replaceable covers

Dimensions(")	Code
780mm L x 450mm H x 220mm Ø	FG6031



FormAlign® P6

The Bolster Roll has many applications from arm supports to trunk support. It's adjustable in size and structural support.

- ✓ Many applications from arm supports to trunk support
- ✓ Adjustable in size and structural support
- ✓ Crumb foam filling effective pressure care
- ✓ Wipeable cover
- ✓ Replaceable covers

Dimensions(")	Code
780mm L x 220mm Ø	FG6053



FormAlign® P1

- ✓ Ideal for support to heels and knees
- ✓ Prevents skin to skin contact
- ✓ Crumb foam filling
- ✓ Moulds to user's body shape
- ✓ Wipeable covers
- ✓ Size: 400mm L x 300mm H

Code: FG6054



FormAlign® P2

- ✓ Ultimately used between the knees
- ✓ Designed to keep the knees from adduction
- ✓ For bed or chair use
- ✓ Wipeable covers
- ✓ Size: 360mm L x 300mm W x 120mm H

Code: FG6055



FormAlign® P3

- ✓ Ideal for supporting upper or lower limbs in elevation
- ✓ Assists with blood circulation
- ✓ Off loads pressure points
- ✓ Wedge cushion, crumb filled
- ✓ Wipeable covers
- ✓ Size: 800mm L x 300mm W x 80/0mm H

Code: FG6056



FormAlign® P5

- ✓ Structural wedge
- ✓ To be used with other "P" crumb foam filled products to add structure to a position
- ✓ Wipeable covers
- ✓ Size: 400mm L x 400mm W x 100/20mm H

Code: FG6058



A photograph showing a person wearing a grey FormAlign harness. Another person's hands are visible, adjusting the black plastic buckle of the harness. The person wearing the harness is wearing a white t-shirt with 'Los Angeles' printed on it. The background is slightly blurred, showing what appears to be an outdoor setting with some equipment.

Belts & Harnesses

Designed with positioning in mind

All the FormAlign belts and harnesses are manufactured to a high standard and are designed with positioning in mind along with comfort and dignity for the user. The range is available in both adult sizes and paediatric and are held in stock for quick delivery.



2 Point Pelvic Belt

Option	Size	Hip Width	Code
Dual Pull Adjustment	Extra Small	5"-9"	FG6005XS
Dual Pull Adjustment	Small	7"-11"	FG6005S
Dual Pull Adjustment	Medium	9"-15"	FG6005M
Dual Pull Adjustment	Large	13"-24"	FG6005L
Centre Pull Adjustment	Extra Small	5"-9"	FG6063XS
Centre Pull Adjustment	Small	7"-11"	FG6063S
Centre Pull Adjustment	Medium	9"-15"	FG6063M
Centre Pull Adjustment	Large	13"-24"	FG6063L
Anti-Tamper	Small	7"-11"	FG6034S
Anti-Tamper	Medium	9"-15"	FG6034M
Anti-Tamper	Large	13"-24"	FG6034L



4 Point Pelvic Belt

Option	Size	Hip Width	Code
Centre Pull Adjustment	Extra Small	5"-9"	FG6062XS
Centre Pull Adjustment	Small	7"-11"	FG6062S
Centre Pull Adjustment	Medium	9"-15"	FG6062M
Centre Pull Adjustment	Large	13"-24"	FG6062L
Anti-Tamper	Small	7"-11"	FG6037S
Anti-Tamper	Medium	9"-15"	FG6037M
Anti-Tamper	Large	13"-24"	FG6037L



Chest Harness with Zip

Size	Shoulder Width	Code
Extra Small	9"-13"	FG6014XS
Small	13"-17"	FG6014S
Medium	17"-20"	FG6014M
Large	21"-24"	FG6014L



Chest Harness without Zip

Size	Shoulder Width	Code
Extra Small	9"-13"	FG6015XS
Small	13"-16"	FG6015S
Medium	17"-20"	FG6015M
Large	21"-24"	FG6015L



Shoulder Harness

Size	Shoulder Width	Code
Extra Small	9"-13"	FG6016XS
Small	13"-17"	FG6016M
Medium	17"-21"	FG6016M
Large	21"-24"	FG6016L



Groin Harness

Size (")	Code
Medium	FG6064M
Large	FG6064L



Calf Straps

Measurements	Code
18" wide x 6" high	FG6040



Ankle Huggers

Size	Code
Small	FG6011S
Medium	FG6011M
Large	FG6011L



D-Rings

Code
FG6060



Universal Strap

Measurements (can be cut)	Code
59 L x 3 W	FG6042

Hardware Tube Clamp

When making life easier is something on the agenda, the hardware clamp product is designed to fit tubular without the need of extra drilling or spacers. The simple yet unique product allows fitting of belts harnesses even easier and hassle free.

- ✔ Can be used in combination with our cam-buckle product
- ✔ Comes complete with bolts
- ✔ Can be purchased as part of a TripSaver fixation kit
- ✔ Seat rails and backrest canes from 17mm-25mm (3/4" to 1")

Code: FG6078



Cam Buckles

Offered in a range of sizes our cam buckles will suit most belt and harness systems on the market. Available in a time saving kit form, it is the perfect solution for the engineer offering the parts readily available for mounting and adjusting the seating system and ensuring there is no compromise and no time lost.

- ✔ 3 widths available
- ✔ Pre-drilled for mounting
- ✔ Available as a trip saver kit solution

Width (")	Code
1	FG6065
1.5	FG6066
2	FG6067



Tri Glides

The FormAlign® Tri Glides allow you to adjust the length of webbing/strap of a belt, making it suit the end user better.

Width (")	Code
1	FG6073
1.5	FG6074
2	FG6075




Harness Retainer Kit


The harness retainer hardware is used to guide the harness over the shoulder and ensure the specific position, prescribed by the clinical team, is maintained for the best result for posture.

- ✔ Maintains harness position
- ✔ Easily adjustable for the prescriber
- ✔ Can be used on third party backrests

Dimensions (")
8 W x 3 H

Code
FGA1101





Notes

Notes

Notes

Notes

Sometimes a little more thought about the life and use of the product can make the world of difference to the user and carer alike.





Contact our team to discuss how FormAlign® can support your organisation:

01473 467 777

sales@formalign.co.uk

formalign.com

Unit 7 Liberator House, Clopton Park, Clopton, Woodbridge, IP13 6QT